

luminfusion



dermalogica PRO



service highlights

LuminFusion is a supercharged skin service which takes luminosity to a whole new level and draws on the use of nanoinfusion with the Dermalogica Pro Pen and Dermalogica's professional-grade exfoliants and masques.

This 45-minute service takes a unique approach with advanced exfoliation stacking, nanoinfusion technology, and LED light therapy.

- Improves skin luminosity.
- Smooths fine lines and wrinkles.
- Instantly reveals radiant skin.

exfoliation stacking

Exfoliation stacking incorporates multiple strengths of exfoliants to target different layers of the skin, providing exceptional results without irritation or down time.

1. **Daily Glycolic Cleanser** unglues and cleanses away the daily buildup of brightness blockers, while maintaining the skin's moisture barrier.
2. **MultiVitamin Power Exfoliant** features retinol, salicylic acid, and lactic acid. This professional strength, lipid-based formula provides faster, more intense exfoliation for optimal smoothness.
3. **Nanoinfusion with the Dermalogica Pro Pen** improves skin texture by enhancing exfoliation and product penetration.
4. **Daily Microfoliant** smooths and brightens with papain, salicylic acid, and rice enzymes.



what is nanoinfusion?

Nanoinfusion is a non-invasive technique utilizing Dermalogica Pro Pen, with stainless steel cone nanotips, which resurfaces the skin and enhances product penetration.

Nanoinfusion helps improve skin texture, boost hydration as well as address signs of premature skin aging. The most immediate, visible change in the skin is a hydrated glow.


the dermalogica technique

The Dermalogica NanoInfusion technique has been designed to amplify service results.

Benefits include:

- Comfortable + Thorough: the tool is moved in a methodical manner.
- Easy To Use: our method allows for fewer stops and starts, with changes in cartridge range and application.
- Health Benefits: encourages lymphatic movement to help the body expel toxins while initiating a positive skin response.
- Even + Complete Treatment: horizontal and vertical passes ensure both.

Depending on your location, there may be restrictions on which licensed professionals are legally authorized to perform certain advanced treatments. Check your local, state, and/or country professional licensing regulations, safety protocols and guidelines to ensure that you are working legally within your scope of practice before performing an advanced treatment.



deep treatment

Light Energy Masque is professional gel masque which is uniquely formulated to amplify the benefits of light-emitting diode technology (LED) combined with red-light therapy.





treatment plans and home care

LuminFusion services can be performed as often as every two weeks and recommended as a series of three for best results.

at-home care

dull skin / dehydration



BioLumin-C Serum
 BioLumin-C Gel Moisturizer
 Invisible Physical Defense SPF 30

hyperpigmentation



PowerBright Dark Spot Serum
 PowerBright Moisturizer SPF 50
 PowerBright Overnight Cream

lines and wrinkles



Dynamic Skin Retinol Serum
 Dynamic Skin Recovery SPF 50
 Phyto Nature Oxygen Cream

pre and post care

Pre-service guidelines are important to follow so you can get the best results from your service.

1. Do not use topical agents that may increase the sensitivity of your skin such as retinoids, exfoliants, topical antibiotics, or, hydroxy acids 5-7 days prior to the service.
2. Avoid IPL/Laser procedures for at least 2 weeks prior to the service.
3. No waxing, depilatory creams, or electrolysis to the area being treated for 5 days prior.
4. Avoid unprotected sun exposure or sunburn at least one week before service.
5. Do not shave the face the day of the procedure to avoid skin irritation.
6. Keep makeup to a minimum day of the service if possible and do not experiment with new products before the service.
7. Avoid alcohol and caffeine 24-48 hours before the service.

Aftercare is critical for achieving your healthiest-looking skin, so we want to make sure you know exactly what to do.

Post Service

1. Try to stay out of excessive sun exposure, such as tanning and sunbathing, for at least 10 days to avoid risk of hyperpigmentation and sensitizing the skin.
2. Do not apply makeup immediately after service. Continue to wear SPF 30+ daily.
3. Before bedtime on the day of service, use a gentle cleanser and tepid water to wash skin.
4. Use serums and the recommended skin care products per your skin therapist's directions.
5. Avoid sweaty exercise, saunas, and swimming for at least 48 hours after service.
6. Avoid immediately using active skin care (including AHA/BHAs, exfoliants, Vitamin A, Retinol, etc.) for at least 3-5 days.

Contraindications

- Active Herpes Simplex
- Active acne
- Skin lesions
- Open cuts, wounds, and abrasions
- Have excessively tanned or sunburned skin from the sun or tanning bed within 14 days
- Seek physician's approval for nanoinfusion if you are pregnant or nursing.

Precautions:

- Wait 14 days after administration of neurotoxin injectables (i.e., Botox).
- Wait 14 days after administration of fillers.
- If using retinol omit use 3-5 days prior to service.
- Avoid alcohol and caffeine 1-2 days prior to service.
- You may experience inflammation, itching, and burning.

dermalogica PRO